



IL4 Antibody, HRP conjugated

Product Code	CSB-PA011659LB01DO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O77762
Immunogen	Recombinant Dog Interleukin-4 protein (25-132AA)
Raised In	Rabbit
Species Reactivity	Dog
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-4 (IL-4) (B-cell stimulatory factor 1) (BSF-1) (Lymphocyte stimulatory factor 1), IL4
Immunogen Species	Dog
Research Area	Immunology
Target Names	IL4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.