



VTG2 Antibody, Biotin conjugated

Product Code	CSB-PA16369D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02845
Immunogen	Recombinant Gallus gallus Vitellogenin-2 protein (26-653AA)
Raised In	Rabbit
Species Reactivity	Gallus gallus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vitellogenin-2 (Major vitellogenin) (Vitellogenin II) [Cleaved into: Lipovitellin-1 (Lipovitellin I) (LVI); Phosvitin (PV); Lipovitellin-2 (Lipovitellin II) (LVII); YGP40], VTG2, VTGII
Immunogen Species	Gallus gallus (Chicken)
Research Area	Others
Target Names	VTG2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.