




MRFAP1L1 Antibody

Product Code	CSB-PA014790EA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96HT8
Immunogen	Recombinant Human MORF4 family-associated protein 1-like 1 protein (1-127AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	MORF4 family-associated protein 1-like 1, MRFAP1L1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MRFAP1L1

Image



Western blot
 All lanes: MRFAP1L1 antibody at 2µg/ml +
 HepG2 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 15 kDa
 Observed band size: 15 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.