



JAK2 Antibody

Product Code	CSB-PA011931LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60674
Immunogen	Recombinant Human Tyrosine-protein kinase JAK2 protein (752-1132AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tyrosine-protein kinase JAK2 (EC 2.7.10.2) (Janus kinase 2) (JAK-2), JAK2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	JAK2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.