





BMP3 Antibody

Product Code	CSB-PA002739LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12645
Immunogen	Recombinant Human Bone morphogenetic protein 3 protein (362-472AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Bone morphogenetic protein 3 (BMP-3) (Bone morphogenetic protein 3A) (BMP-3A) (Osteogenin), BMP3, BMP3A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	BMP3

Image



Western blot
All lanes: BMP3 antibody at 2µg/ml
Lane 1: SW1990 whole cell lysate
Lane 2: 293T whole cell lysate
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 54 kDa
Observed band size: 54 kDa



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA002739LA01HU at dilution of 1:50