




PGC Antibody

Product Code	CSB-PA11299A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20142
Immunogen	Recombinant Human Gastricsin protein (17-388AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gastricsin (EC 3.4.23.3) (Pepsinogen C), PGC
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PGC

Image



Western Blot

Positive WB detected in: Mouse stomach tissue
 All lanes: PGC antibody at 3.4µg/ml
 Secondary
 Goat polyclonal to rabbit IgG at 1/50000 dilution
 Predicted band size: 43, 35 kDa
 Observed band size: 43, 35 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.