





MKI67 Antibody

Product Code	CSB-PA11617A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46013
Immunogen	Recombinant Human Proliferation marker protein Ki-67 protein (3120-3256AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proliferation marker protein Ki-67 (Antigen identified by monoclonal Ki-67) (Antigen KI-67) (Antigen Ki67), MKI67
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MKI67

Image



Immunohistochemistry of paraffin-embedded human endometrial cancer using CSB-PA11617A0Rb at dilution of 1:100



Immunohistochemistry of paraffin-embedded human bladder cancer using CSB-PA11617A0Rb at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.