



rpiA Antibody, HRP conjugated

Product Code	CSB-PA540419DB01ENT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A7Z0
Immunogen	Recombinant Escherichia coli Ribose-5-phosphate isomerase A protein (1-219AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribose-5-phosphate isomerase A (EC 5.3.1.6) (Phosphoriboisomerase A) (PRI), rpiA, ygfC
Immunogen Species	Escherichia coli (strain K12)
Research Area	Signal Transduction
Target Names	rpiA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.