



esxB Antibody, Biotin conjugated

Product Code	CSB-PA363658YD01MVZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P9WVK4
Immunogen	Recombinant Mycobacterium tuberculosis ESAT-6-like protein EsxB protein (2-100AA)
Raised In	Rabbit
Species Reactivity	Mycobacterium tuberculosis
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ESAT-6-like protein EsxB (10 kDa culture filtrate antigen CFP-10) (Secreted antigenic protein MTSA-10), esxB, cfp10 lhp mtsA10
Immunogen Species	Mycobacterium tuberculosis (strain CDC 1551 / Oshkosh)
Research Area	Others
Target Names	esxB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.