




TUBB Antibody

Product Code	CSB-PA025318OJ01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07437
Immunogen	Peptide sequence from Human Tubulin beta chain protein
Raised In	Rabbit
Species Reactivity	Human, Zebrafish
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.1mol/L NaCl, 15mmol/L NaN ₃ , PH7.2
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tubulin beta chain (Tubulin beta-5 chain), TUBB, TUBB5
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	TUBB

Image



Western blot

All lanes: TUBB antibody at 6µg/ml

Lane 1: PC-3 whole cell lysate

Lane 2: 293T whole cell lysate

Lane 3: A549 whole cell lysate

Lane 4: HepG2 whole cell lysate

Lane 5: SH-SY5Y whole cell lysate

Lane 6: Hela whole cell lysate

Lane 7: Zebrafish lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 50 kDa

Observed band size: 50 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.