



Tpsb2 Antibody, HRP conjugated

Product Code	CSB-PA024128LB01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21845
Immunogen	Recombinant Mouse Tryptase beta-2 protein (32-276AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tryptase beta-2 (Tryptase-2) (EC 3.4.21.59) (Mast cell protease 6) (mMCP-6), Tpsb2, Mcpt6
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	Tpsb2