




ROR1 Antibody

Product Code	CSB-PA020067LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01973
Immunogen	Recombinant Human Inactive tyrosine-protein kinase transmembrane receptor ROR1 protein (30-391AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:10-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inactive tyrosine-protein kinase transmembrane receptor ROR1 (Neurotrophic tyrosine kinase, receptor-related 1), ROR1, NTRKR1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ROR1

Image



Immunohistochemistry of paraffin-embedded human prostate tissue using CSB-PA020067LA01HU at dilution of 1:10

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.