



Ces1c Antibody, FITC conjugated

| | |
|----------------------------|--------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-PA338557YC01MO |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P23953 |
| Immunogen | Recombinant Mouse Carboxylesterase 1C protein (19-550AA) |
| Raised In | Rabbit |
| Species Reactivity | Mouse |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Carboxylesterase 1C (EC 3.1.1.1) (Liver carboxylesterase N) (Lung surfactant convertase) (PES-N), Ces1c, Es1 |
| Immunogen Species | Mus musculus (Mouse) |
| Research Area | Others |
| Target Names | Ces1c |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |