



Ces1c Antibody, FITC conjugated

Product Code	CSB-PA338557YC01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23953
Immunogen	Recombinant Mouse Carboxylesterase 1C protein (19-550AA)
Raised In	Rabbit
Species Reactivity	Mouse
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Carboxylesterase 1C (EC 3.1.1.1) (Liver carboxylesterase N) (Lung surfactant convertase) (PES-N), Ces1c, Es1
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Ces1c