



# Ces1c Antibody, FITC conjugated

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA338557YC01MO   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  |
| <b>Uniprot No.</b>         | P23953   |
| <b>Immunogen</b>           | Recombinant Mouse Carboxylesterase 1C protein (19-550AA)   |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Mouse  |
| <b>Form</b>                | Liquid   |
| <b>Conjugate</b>           | FITC   |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4                             |
| <b>Purification Method</b> | >95%, Protein G purified   |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Alias</b>               | Carboxylesterase 1C (EC 3.1.1.1) (Liver carboxylesterase N) (Lung surfactant convertase) (PES-N), Ces1c, Es1 |
| <b>Immunogen Species</b>   | Mus musculus (Mouse)   |
| <b>Research Area</b>       | Others   |
| <b>Target Names</b>        | Ces1c  |