



Recombinant Schizosaccharomyces pombe Zinc finger protein zpr1 (zpr1)

Product Code	CSB-BP026617SXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O13724
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MAEEKKEELF TSIGNAAQNV STAEDREGNG VQEVESLCME CGKNGTTKLL LTVIPIYFREV VLMSFECPHC GFKNAQVQHA ETIQPEGSKI TFHVEDKEDL NRTVVKSQEA IVSIPEIQLE IPGRLGQLTT IEGILSNVVD DLSKEQESRK ESAPQLYDQI NAFIEKVNSL RSGSVPFTIT VDDITGNSWI EMKPGRDGDR WSQVSYKRTL EQNTKLGLVD TDQPEDVKTQ TNNASNTLKH DATAVEVDPN EVHTFHATCP SCSHQCDTHM KLLDIPHFKE VIIMSTVCDR CGYRSNEVKT GGEIPPKGRK ITLKVMDAED LSRDILKSET ASLKIPELGL DLFPGTLGGR FTTIEGLLAQ VYDELYGRVF SQETDSMTPE QVANWQQFLC NLTAAREGAT QFTLILDDPL SQSYLQNYA PDPDPNMTIE EYERSFQVNE ELGLNDMKTE NYEKDGGKK
Source	Baculovirus
Target Names	zpr1
Protein Names	Recommended name: Zinc finger protein zpr1
Expression Region	1-459
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.