




Lox Antibody

Product Code	CSB-PA013038LA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16636
Immunogen	Recombinant Rat Protein-lysine 6-oxidase protein (163-411AA)
Raised In	Rabbit
Species Reactivity	Rat, Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein-lysine 6-oxidase (EC 1.4.3.13) (Lysyl oxidase), Lox
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Signal Transduction
Target Names	Lox

Image



Western blot

All lanes: Lox antibody at 2µg/ml + K562 whole cell lysate
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 47 kDa
Observed band size: 47 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.