



# Rabbit anti-Chicken Yolk Immunoglobulin Antibody; Biotin conjugated

<b>Product Code</b>	CSB-PA00410H0Rb
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Immunogen</b>	Chicken Yolk Immunoglobulin
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Chicken
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	Caprylic Acid Ammonium Sulfate Precipitation purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Product Type</b>	Secondary Antibody
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.