



Rabbit anti-Chicken Yolk Immunoglobulin Antibody;HRP conjugated

Product Code	CSB-PA00410F0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Chicken Yolk Immunoglobulin
Raised In	Rabbit
Species Reactivity	Chicken
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Caprylic Acid Ammonium Sulfate Precipitation purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Secondary Antibody
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.