



Rabbit anti-Chicken Yolk Immunoglobulin Antibody;HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA00410F0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Immunogen | Chicken Yolk Immunoglobulin |
| Raised In | Rabbit |
| Species Reactivity | Chicken |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | Caprylic Acid Ammonium Sulfate Precipitation purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Product Type | Secondary Antibody |