





NT5C1A Antibody

Product Code	CSB-PA866304YA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXI3
Immunogen	Recombinant Human Cytosolic 5'-nucleotidase 1A protein (1-368AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:10-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytosolic 5'-nucleotidase 1A (cN1A) (EC 3.1.3.5) (Cytosolic 5'-nucleotidase IA) (cN-I) (cN-IA), NT5C1A
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	NT5C1A

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA866304YA01HU at dilution of 1:10



Immunofluorescent analysis of A549 cells using CSB-PA866304YA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.