



Recombinant Arabidopsis thaliana Probable ribosome biogenesis protein RLP24 (At2g44860)

Product Code	CSB-BP020543DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O22165
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MRLEKCFWFCSTIYPGHGIQFVRNDAKIFR FCRSKCHKNF KMKRNPRKVK WTKAFRAAHG KDMTKDTTFE FEKKNRPER YDRNVTENTL MAIKKI AKIR TAREAKHIEN RLKPNKQKKL NDDMKELDQN IHMVQAPGAQ KIKVDVSEKK SVQNEAMEE
Source	Baculovirus
Target Names	At2g44860
Protein Names	Recommended name: Probable ribosome biogenesis protein RLP24
Expression Region	1-159
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.