



Recombinant Zea mays S-adenosylmethionine decarboxylase proenzyme (SAMDC)

| | |
|--------------------------|---|
| Product Code | CSB-EP001658ZAX-B |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | O24575 |
| Product Type | Recombinant Protein |
| Immunogen Species | Zea mays (Maize) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MAVLSAADAS PVSAIGFEGY EKRLEITFSE APVFVDPHGR GLRALSRAQI DSVLDLARCT IVSELSNKDF DSYVLSE |
| Source | E.coli |
| Target Names | SAMDC |
| Protein Names | Recommended name: S-adenosylmethionine decarboxylase proenzyme Short name= AdoMetDC Short name= SAMDC EC= 4.1.1.50 Cleaved into the following 2 chains: 1. S-adenosylmethionine decarboxylase alpha chain 2. S-adenosylmethion |
| Expression Region | 1-77 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |