



Recombinant Methanothermobacter thermautotrophicus Signal recognition particle 19 kDa protein (srp19)

Product Code	CSB-EP022674MSR
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O26267
Product Type	Recombinant Protein
Immunogen Species	Methanothermobacter thermautotrophicus (strain ATCC 29096 / DSM 1053 / JCM 10044 / NBRC 100330 / Delta H) (Methanobacterium thermoautotrophicum)
Purity	>85% (SDS-PAGE)
Sequence	MNLIIWPTYL DSRKSRSEGR RVPLEYAVES PTASEILRAA RKLQLEASME SDRAYPPSWW ESSGRVVVEY NGKKSELLPK IARLVRSSRK R
Source	E.coli
Target Names	srp19
Protein Names	Recommended name: Signal recognition particle 19 kDa protein Short name= SRP19
Expression Region	1-91
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.