



Recombinant *Caenorhabditis elegans* Acid ceramidase (asah-1)

| | |
|--------------------------|--|
| Product Code | CSB-YP002169CXY |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | O45686 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Caenorhabditis elegans</i> |
| Purity | >85% (SDS-PAGE) |
| Sequence | KHVE LPAPFKDHCI LDDKQONLYDP SKQFDIKWYD VNLDLPPSER WVQIATANKE HIADLIGVLI NLITPWFNPA IDFVDDVFGD LAPKLAQPYR DEIFSIANAT GIPLGQITMY NIFYEIFTVC TSVIAQDKDG HVFHARNLDF GLFMGWDPVL HDWQISQKLR KMIINVNWLK DGKLLYKSNN FAGYIGIYNG LKPNAFSLTA DDRFQLVGGY YGILKWVFGLEADGKWMSWL ARETLETKTT YLDAKEHLMN TPMLSPVYFI LGGSKKDEGC IARS�DKTA LLEMATSPH GWYLLETNYD QGTEDLYLDD RDTPGFRCMD KLTQKNVGFE GIFNVLSRT NLNKLTTYTV LMSVETSRFE TILQSCPGECPYW |
| Source | Yeast |
| Target Names | asah-1 |
| Protein Names | Recommended name: Acid ceramidase EC= 3.5.1.23Alternative name(s): Acylsphingosine amidohydrolase 1 |
| Expression Region | 17-393 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |