



Recombinant 60S ribosomal protein L18 (rpl-18)

Product Code	CSB-MP020154CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O45946
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MGIDINHKHD RVARRTAPKS ENPYLRLLSK LYAFLARRTG EKFNAIVLKR LRMSRRNRQP LSLAKLARAV QKAGNENKTV VTLSTVTDDA RLYTVPKISV AALHVTEGAR ARILAAGGEI ITLDQLALKS PKGENTVFLQ GPRSAREAEK HFGPAPGVPH SHTKPYVRSK GRKFERARGR RASRAYKN
Source	Mammalian cell
Target Names	rpl-18
Protein Names	Recommended name: 60S ribosomal protein L18
Expression Region	1-188
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.