



Recombinant *Drosophila melanogaster* Flotillin-2 (Flo-2)

Product Code	CSB-EP008728DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O61492
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MGNIHTTGPN EALIVSGGCC GSTKKRTIVG GWAWAWWLVT DVQRSLNVM TLNPMCENVE TSQGVPLTVT GVAQCKIMKS SSKYQTDYHN DEADELLGTA SEQFLGKSVK EIKQTLQTL EGHRLAILGT LTVVEVYKDR DQFAALVREV AAPDVGRMGI EILSFTIKDV YDDVQYLA SL GKAQTAVVKR DADAGVAEAN RDAGIREAEC EKSAMDVKYS TDTKIEDNTR MYKLQKANFD QEINTAKAES QLAYELQAAK IRQRIRNEEI QIEVVERRKQ IEIESQEVQR KDRELTGTVK LPAEAEAFRL QTLAQAKQCQ TIEGARAEAE RIRKIGSAEA HAIELVGKAE AERM RMKAHV YKQYGDAAIM NIVLESLPKI AAEVAAPLAK TDEIVLIGGN DNITNDVTRL VAQLPPSINA LTGVDLSKVL SKIPGAKA
Source	E.coli
Target Names	Flo2
Protein Names	Recommended name: Flotillin-2
Expression Region	1-438
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.