



Recombinant *Oryza sativa* subsp. japonica Proteasome subunit beta type-1 (PBF1)

Product Code	CSB-YP018876OFG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O64464
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MSRRGDWVYE NNGGTCVAIA GADYCVVAAD TRLSVGYNIL TRDHSKICEL ADKCALASSG FQGDIKALHK NLAARELLYQ HQHNKRMSCP AMAQLLSNTL YYKRFFPYA FNVLGGLDSE GKGCVFYTYDA VGSYERTGYS AQGTGSALIM PVLNQLKSP SPLLLPARDA VTPLSETEAV DLVKDFASA TERDIYTGDK LEIVVINKAG TKREYIDLRK D
Source	Yeast
Target Names	PBF1
Protein Names	Recommended name: Proteasome subunit beta type-1 EC= 3.4.25.1 Alternative name(s): 20S proteasome alpha subunit F 20S proteasome subunit beta-6
Expression Region	1-221
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.