



Recombinant *Xenopus laevis* Guanine nucleotide-binding protein subunit alpha-14 (gna14)

Product Code	CSB-BP009586XBE
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O73819
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	>85% (SDS-PAGE)
Sequence	MAGCCLSAEE KESQRINAEI EKQLRRDKRD ARRELKLLLL GTGESGKSTF IKQMRIIHGS GYTDEDRKGF TKLVYQNIPT SMQSMIRAMD TLRIQYTSEQ NMENALVIRE VEVDKVSSLE RKHVEAIKKL WEDEGIQECY DRRREYQLSD STKYYLSDID RISNPGFIPT QQDVLRVVVP TTGIIYFPD LENIIFRMVD VGGQRSERRK WIHCFENVTS IIFLVALSEY DQVLAECDNE NRMEESKALF KTIITYPWFQ NSSVILFLNK KDLLQEKIMY SHLIDYFPEF TGPKQDSQAA RDFILKLYQD QNPDKEKVIY SHFTCATDTE NIRFVFAAVK DTILQLNLRE FNLV
Source	Baculovirus
Target Names	gna14
Protein Names	Recommended name: Guanine nucleotide-binding protein subunit alpha-14 Short name= G alpha-14 Short name= G-protein subunit alpha-14
Expression Region	1-354
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.