



Recombinant Human Septin-4 (SEPTIN4)

Product Code	CSB-MP021027HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O43236
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MDRSLGWQGN SVPEDRTEAG IKRFLEDTTD DGELSKFVKD FSGNASCHPP EAKTWARPQ VPEPRPQAPD LYDDDLEFRP PSRPQSSDNQ QYFCAPAPLS PSARPRSPWG KLDPYDSSSED DKEYVGFATL PNQVHRKSVK KGFDFTLMVA GESGLGKSTL VNSLFLTDLY RDRKLLGAE E RIMQTVEITK HAVDIEEKGV RLRLTIVDTP GFGDAVNTE CWKPVAEYID QQFEQYFRDE SGLNRKNIQD NRVHCCLYFI SPFGHGLRPL DVEFMKALHQ RVNIVPILAK ADTLTPPEVD HKKRKIREEI EHFGIKIYQF PDCDSDEDED FKLQDQALKE SIPFAVIGSN TVVEARGRRV RGRLYPWGIV EVENPGHCDF VKLRTMLVRT HMQDLKDVTR ETHYENYRAQ CIQSMTRLVV KERNRNKLTR ESGTDFPIPA VPPGTDPETE KLIREKDEEL RRMQEMLHKI QKQMKENY
Source	Mammalian cell
Target Names	SEPT4
Protein Names	Recommended name: Septin-4 Alternative name(s): Apoptosis-related protein in the TGF-beta signaling pathway Short name= ARTS Bradeion beta Brain protein H5 CE5B3 beta Cell division control-related protein 2 Short name=
Expression Region	1-478
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.