



Recombinant *Oryza sativa* subsp. japonica Tubulin gamma-2 chain (TUBG2)

Product Code	CSB-BP025336OFG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O49068
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MPREIITIQV GQCGNQIGME FWKQLCLEHG IGKDGLLEDF ATQGGDRKDV FFYQADDQHY IPRALLVDLE PRVINGIQNS EYRNLYNHEN IFVAEHGGGA GNNWASGYHQ GEQVVDDIMD MIDREADGSD SLEGFVLCHS IAGGTGSGMG SYLLETLNDR YSKKLVQTY S VFPNQMETSD VVVQPYNLL TLKRLTLNAD CVVVDNTAL NRIAVERLHL ANPTFAQTNS LVSTVMSAST TTLRYPGYMN NDLVGLLASL IPTPRCHFLM TGYTPLTVER QVNMIRKTTV LDVMRRLQT KNIMVSSYAR NKEASQAKYI SILNIIQGEV DPTQVHESLQ RIRERKLVNF IEWGPASIQV ALSRKSPYVQ TTHRVSGLML ANHTSIRHLF SKCLGQYEKL RKKQAFLDNY RKFPMFEDND LSEFDESREI IESLVDEYKA CESPDYIKWG MEDAGEANVA AALDSKLVV
Source	Baculovirus
Target Names	TUBG2
Protein Names	Recommended name: Tubulin gamma-2 chain Alternative name(s): Gamma-2-tubulin
Expression Region	1-469
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.