




Recombinant Oncorhynchus mykiss Transforming growth factor beta-1 (TGFB1)

Product Code	CSB-YP023446OEI
Relevance	Likely to be an important cytokine regulating immune response. May also have a role in other physiological systems.
Abbreviation	Recombinant Oncorhynchus mykiss TGFB1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O93449
Alias	Short name: TGF-beta-1
Product Type	Recombinant Protein
Immunogen Species	Oncorhynchus mykiss (Rainbow trout) (Salmo gairdneri)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	QTTTEEICSDKSESCCVRKLYIDFRKDLGWKWIHEPTGYFANYCIGPCTYIWNT ENKYSQVLALYKHHNPGASAPCCVPQVLEPLPIIYYVGRQHKVEQLSNMIVK SCRCS
Research Area	Signal Transduction
Source	Yeast
Target Names	tgfb1
Protein Names	Recommended name: Transforming growth factor beta-1 Short name= TGF-beta-1
Expression Region	271-382aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	15.0kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.