



# Recombinant Coturnix coturnix japonica Estrogen receptor beta (ESR2)

<b>Product Code</b>	CSB-MP007831DXJ
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O93511
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Coturnix japonica (Japanese quail) (Coturnix coturnix japonica)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAFCSPAMMN YNIASNFGDS ESASVRQTSS PSELLWSAPGH LSPLTLHCQL SLLYAEQPKS PWCEARPLEP VLPVSRETLK RKTNGSDCTS PIASNPGSKR DAHFCVCS D YASGYHYGVW SCEGCKAFFK RSIQGHNDYI CPATNQCTID KNRRKSCQAC RLRKCYEVGM MKCGSRRERC GYRILRRHRN SEDCMGKTKK YNEAATRVKE ILLSTVSPEQ FVLTLLAEAP PNVLVS RPSK PFTEASMMMS LTKLADKELV HMIGWAKKIP GFIDLSLYDQ VRLLESCWLE VLMIGLMWRS IDHPGKLIFA PDLVLD RDEG KCVEGILEIF DMLLAMTSRF RELKLQHKEY LCVKAMILLN SSMFPLSAEE PESNRKLHHL LNVVTDALVW VIAKSGIPSQ QQTTRLANLL MLLSHVRHAS NKGMEHLLSM KCKNVVPVYD LLEMLNAHT LRGQRKSPVT HPDFEQVSHF QV
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ESR2
<b>Protein Names</b>	Recommended name: Estrogen receptor beta Short name= ER-beta Alternative name(s): Nuclear receptor subfamily 3 group A member 2
<b>Expression Region</b>	1-472
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.