



# Recombinant Xiphophorus helleri Cellular tumor antigen p53 (tp53)

<b>Product Code</b>	CSB-BP024077XBJ
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O57538
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Xiphophorus helleri (Green swordtail)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEEADLTPL SQDTFHDLWN NVFLSTENES LAPPEGLLSQ NMDFWEDPET MQETKNVPTA PTVPAISNYA GEHGFNLEFN DSGTAKSVTS TYSVKLGKLF CQLAKTTPIG VLVKEEPPQG AVIRATSVYK KTEHVGVEVK RCPHHQSEDL SDNKSHLIRV EGSQLAQYFE DPNTRRHSVT VPYERPQLGS EMTTILLSFM CNSSCMGGMN RRPILTILT ETTEGEVLGR RCFEVRVCAC PGRDRKTEEG NLEKSGTKQT KKRKSAPAPD TSTAKKSKSA SSGEDEDKEI YTLSIRGRNR YLWFKSLNDG LELMDKTGPK IKQEIPAPSS GKRLKGGSD SD
<b>Source</b>	Baculovirus
<b>Target Names</b>	tp53
<b>Protein Names</b>	Recommended name: Cellular tumor antigen p53 Alternative name(s): Tumor suppressor p53
<b>Expression Region</b>	1-342
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.