



# Recombinant Schizosaccharomyces pombe Nucleolar protein 56 (nop56)

<b>Product Code</b>	CSB-YP015939SXV
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O94514
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MADYLLYESA TGYSLFDVVG ADQIAAKTKE VQLSLQDISK FGKVVQLRSF IPFKNAHAL ENANDISEGV LNDLKNFLE LNLPKASKKK KVS LGVQDKN LATSİKSEİD AİECDTSELT QDLLRGİRFH GDKLLKQLSP GDFERAQLGL GHSYSRAKVK FNVNRNDNMI IQAİAILDQL DKDİNTFAMR MKEWYSWHFP ELSKİVGDNY KYAVİVTLVG DKTTİNDEML HDLAAVVDDD KDİAQSIİNA GKVSİMGQDIS EİDLENİLSF AERVİKLSNY RKQLHNYLVQ KMNVVAPNLA ELİGEMVGAR LİSHAGSLTN LSKCPASTVQ İLGAEKALFR ALKTRGNTPK YGİİYHSSFI GKAGAKNKGR İSRFLANKCS İASRIDNFSD APTTAFGQVL RRQVEERLNF FDTGVAPTRN SİAMAEAYEK ALSSVNİDGD EEVDİDVEET VETİSEKPSK KEKKDKKEKK KEKSİKKKRSİ DDASEEVKKS KKKKKSH
<b>Source</b>	Yeast
<b>Target Names</b>	nop56
<b>Protein Names</b>	Recommended name: Nucleolar protein 56 Alternative name(s): Ribosome biosynthesis protein sik1
<b>Expression Region</b>	1-497
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.