



Recombinant Spalax ehrenbergi Somatotropin (GH1)

Product Code	CSB-EP009407FJV-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O70615
Product Type	Recombinant Protein
Immunogen Species	Spalax ehrenbergi (Middle East blind mole rat) (Nannospalax ehrenbergi)
Purity	≥85% (SDS-PAGE)
Sequence	FPAM PLSNLFANAV LRAQHLHQLA ADTYKEFERA YIPEGQRYSI QNAQA AFCFS ETIPAPTGKE EAQQRSDMEL LRFSLLLIQS WLGPVQFLSR VFTNSLVFGT SDRVFEKLD LEEGIQALMR ELEDGSLRAG QLLKQTYDKF DTNMRSDDAL LKNYGLLS CF KKD LHKAEY LRVMKCRRFV ESSCAF
Source	E.coli
Target Names	GH1
Protein Names	Recommended name: Somatotropin Alternative name(s): Growth hormone
Expression Region	27-216
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.