



Recombinant Schizosaccharomyces pombe Transcription initiation factor TFIID subunit 6 (taf6)

Product Code	CSB-EP023095SXV-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O74462
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSLTVWNIES IKDVAEMLGI GNLADEPAAA IAMDLEYRIH QVVQEATKFM VHSKRTVLTS ADISSALRTL NVEPLYGFNN SRPLEFHEAA VGAGQNSLYY LDDEEVDFEK IINAPLPKVP RNISYSAHWL AIEGVQPAIP QNPTPSDHTV GEWASKGTSG VMPGASTAAK EARNGVTSM D NVEIKPLVRH VLSKELQLYF ERITSALLDE TNVELRDAAL SSLRDDPGLH QLLPYFIMFL SDSVTRNLGN LVVLTTLMHM AWALLDNP NL FVEPYVQQLM PSILTCLVAK RLGSDPNNHE HYALRDAAF LLGIVCDRFG NVYYTLKPRV TRTALKAF LD NTKPYSTHYG AIKGLKTMGK EAIRVLVVPN IKVYEVLVRK TLEKGNEEEI YEANKCMDAL YDALLLRDD QLPNQRTLPP NASGLLEKNV GSLMAEKIMK ENDTSLLLGL LE
Source	E.coli
Target Names	taf6
Protein Names	Recommended name: Transcription initiation factor TFIID subunit 6 Alternative name(s): TBP-associated factor 50 kDa Short name= TAFII-50 Short name= TAFII50 TBP-associated factor 6
Expression Region	1-452
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.