



Recombinant Human Estrogen-related receptor gamma (ESRRG)

| | |
|--------------------------|---|
| Product Code | CSB-EP007837HU |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P62508 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MDSVELCLPE SFSLHYEEL LCRMSNKDRH IDSSCSSFIK TEPSSPASLT DSVNHHSPPG SSDASGSYSS TMMNGHQGLD SPPLYPSAPI LGGSGPVRKL YDDCSSTIVE DPQTKCEYML NSMPKRLCLV CGDIASGYHY GVASCEACKA FFKRTIQGNI EYSCPATNEC EITKRRRKSC QACRFMKCLK VGMLKEGVRL DRVRGGRQKY KRRIDAENSP YLNPQLVQPA KKPYNKIVSH LLVAEPEKIY AMPDPTVPDS DIKALTTLCD LADRELVVII GWAKHIPGFS TSLADQMSL LQSAWMEILI LGVVYRSLSF EDELVYADDY IMDEDQSKLA GLLDLNNAIL QLVKKYKSMK LEKEEFVTLK AIALANSDSM HIEDVEAVQK LQDVLHEALQ DYEAGQHMED PRRAGKMLMT LPLLRTSTK AVQHFYNIKL EGKVPMHKLF LEMLEAKV |
| Source | E.coli |
| Target Names | ESRRG |
| Protein Names | Recommended name: Estrogen-related receptor gamma Alternative name(s): ERR gamma-2 Estrogen receptor-related protein 3 Nuclear receptor subfamily 3 group B member 3 |
| Expression Region | 1-458 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |