



# Recombinant Human Cellular nucleic acid-binding protein (CNBP)

<b>Product Code</b>	CSB-YP005637HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P62633
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	SSNECFKCG RSGHWARECP TGGGRGRGMR SRGRGGFTSD RGFQFVSSSL PDICYRCGES GHLAKDCDLQ EDACYNCGRG GHIAKDCKEP KREREQCCYN CGKPGHLARD CDHADEQKCY SCGEFGHIQK DCTKVKCYRC GETGHVAINC SKTSEVNCYR CGESGHLARE CTIEATA
<b>Source</b>	Yeast
<b>Target Names</b>	CNBP
<b>Protein Names</b>	Recommended name: Cellular nucleic acid-binding protein Short name= CNBP Alternative name(s): Zinc finger protein 9
<b>Expression Region</b>	2-177
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene encodes a nucleic-acid binding protein with seven zinc-finger domains. The protein has a preference for binding single stranded DNA and RNA. The protein functions in cap-independent translation of ornithine decarboxylase mRNA, and may also function in sterol-mediated transcriptional regulation. A CCTG expansion in the first intron of this gene results in myotonic dystrophy type 2. Multiple transcript variants encoding different isoforms have been found for this gene.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.