



Recombinant Human Small EDRK-rich factor 2 (SERF2)

Product Code	CSB-MP021042HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P84101
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MTRGNQRELA RQKNMKKQSD SVK GKRRDDG LSAAARKQRD SEIMQQKQKK ANEKKEEPK
Source	Mammalian cell
Target Names	SERF2
Protein Names	Recommended name: Small EDRK-rich factor 2 Alternative name(s): Gastric cancer-related protein VRG107 Protein 4F5-related Short name= 4F5rel Short name= h4F5rel
Expression Region	1-59
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.