



Recombinant Arabidopsis thaliana Cytoplasmic tRNA 2-thiolation protein 1 (NCS6)

Product Code	CSB-MP002426DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O64862
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MEAKNKKAVA SRLCCLCNLR RPVLKRPKTL QQICRECFYE VFEEIHQVI VQNRLFKSGE RVAIGASGGK DSTVLAYVLS ELNRRHNYGL DFLLSIDEG ITGYRDDSL E TVKRNEVQYG LPLKIVSYKD LYGWTMDEIV KMIGLKNCT FCGVFRRQAL DRGAALLKVE KLVTGHNADD IAETVLLNIL RGDARLSRC TSITTGEDGP IPRCKPFKYT YEKEIVMYAY FKKLDYFSTE CIYSPNAYRG FAREFIKDLE RIRPRAILDI IKSGEDFRIA TTTKMPEQGT CERCYISSQ KWCKACVLE GLNRGLPKMG IGRPRGVNGD HNKETKKPGS VAKSIESKQC GSLDF
Source	Mammalian cell
Target Names	NCS6
Protein Names	Recommended name: Cytoplasmic tRNA 2-thiolation protein 1 EC= 2.7.7.- Alternative name(s): Cytoplasmic tRNA adenylyltransferase 1
Expression Region	1-355
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.