



Recombinant Schizosaccharomyces pombe DNA primase large subunit (pri2)

Product Code	CSB-BP018681SXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O74761
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MFRTTKSRVV EKNQNFSSAS FESLSYPTRL NFKKPPVSE ISIEEFETWA IDRLVVLGEI ESAMLRNKT NELNGIHKI LDKRLPMSSN MARTVKGNL DEERRKDHYS HFILRLAFSR SDELRI RFLR AESTLFRFRF VNEETRERQA FIDSLDFQWE SVSQEEKDVF YEKLQASSQR NIENESFFKV PFFKVPDLVE RRAV FVHKGY AYVPFSEQVS LVVEEFEENL KKALESTAKS LPRLEDDRL LPILNHL SKG FVAPESSIVA PKSGAITAGQ VDSFVSHFPL CAKHLHEVLK RDKHLRHFGR LQYGLFLKDI GLSVDEALVF WRKSFTNVTE DKFNKEYRYN IRHTYGLEGN RKNYKGYNCQ QILTGPQLGP GDAHGC PYRT YSVNNLVSAL ASMG IAPDSA GCREVVEAVK ARHYHI ACTR VFELTHPNVG SLEESIHPL QYFQLSLES
Source	Baculovirus
Target Names	pri2
Protein Names	Recommended name: DNA primase large subunit EC= 2.7.7.-
Expression Region	1-459
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.