



Recombinant Ribosomal RNA-processing protein 8 (T07A9.8)

Product Code	CSB-EP020531CXY-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O44410
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MGKKRKITDE KDAQHVPAEK REKVENWLKK STEKPTSSQS DAEKKKKRPW RNKVRKLAAK KAAADKksen PEEPPLILEP KSSSDENTKK KRKRGPKkkk FKPEVAGKAA ETENDDVAAA PEEADPIAEA KKRLDAGRFR FLNEKLYTCT GSEAFDFFKE DPTAFDLYHK GFADQVKKWP NHPLREIIRW LQSKPDQQSV FDLGCGEAKI AEAUGEKHKI RSFDLVAVND RVESCDMSKL PAEDSSADIV IYCLSLMGTN LYDFIREARR VLKIGGILKI AEVTSRFVSI KQFCEAITKM GFEQSHRREL TDYFMMMEFK KVEKVEQKRP YGLKLPCLY KKR
Source	E.coli
Target Names	T07A9.8
Protein Names	Recommended name: Ribosomal RNA-processing protein 8 EC= 2.1.1.-
Expression Region	1-343
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.