



# Recombinant Human Ribonuclease P protein subunit p40 (RPP40)

<b>Product Code</b>	CSB-BP020351HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O75818
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MATLRRLREA PRHLLVCEKS NFGNHKSRHR HLVQTHYYNY RVSFLIPECG ILSEELKNLV MNTGPYYFVK NLPLHELITP EFISTFIKKG SCYALTYNTH IDEDNTVALL PNGKLILSLD KDTYEETGLQ GHPSQFSGRK IMKFIVSIDL MELSLNLDK KYERISWSFK EKKPLKFDL LAWHTGSEE STMMSYFSKY QIQEHQPKVA LSTLRDLQCP VLQSSELEGT PEVSCRALEL FDWLGAVFSN VDLNNPEPNF ISTDCCPEPS TVVAKAYLCT ITGFILPEKI CLLLEHLCHY FDEPKLAPWV TLSVQGFADS PVSWEKNEHG FRKGGEHLYN FVIFNNQDYW LQMAVGANDH CPP
<b>Source</b>	Baculovirus
<b>Target Names</b>	RPP40
<b>Protein Names</b>	Recommended name: Ribonuclease P protein subunit p40 Short name= RNaseP protein p40 EC= 3.1.26.5 Alternative name(s): RNase P subunit 1
<b>Expression Region</b>	1-363
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.