



Recombinant *Coryphopterus nicholsi* L-lactate dehydrogenase A chain (Idha)

Product Code	CSB-MP012832DXE
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O93401
Product Type	Recombinant Protein
Immunogen Species	Rhinogobiops nicholsii (Blackeye goby) (<i>Coryphopterus nicholsi</i>)
Purity	>85% (SDS-PAGE)
Sequence	STKEKLISH VSKEEAVGSR NKVTVVGVGM VGMASAISIL LKDLCDLAL VDVMEDKCLKG EVMDLQHGSL FLKTHKIVGD KDYSVTANSR VVVVVTAGARQ QEGESRLNLV QRNVNIFKFI IPNIVKYSPN CILMVVSNPV DILTYVAWKL SGFPRHRVIG SGTNLD SARF RHIMGEKHL HPSSCHGWIV GEHGDSSVPV WSGVNVAGVS LQTLNPKMGA EGDSENWKAV HKMVVDGAYE VIKCLKGYTSW AIGMSVADLV ESIVKNLHKV HPVSTLVKGM HGVKDEVFLS VPCVLGNSGL TDVIHMTLKA DEEKQLVKS A ETLWGVQKEL TL
Source	Mammalian cell
Target Names	Idha
Protein Names	Recommended name: L-lactate dehydrogenase A chain Short name= LDH-A EC= 1.1.1.27
Expression Region	2-332
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.