



# Recombinant Schizosaccharomyces pombe Pre-mRNA-splicing factor 18 (prp18)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-MP018760SXV   |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.   |
| <b>Uniprot No.</b>       | O94406  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)   |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MDFLKEEIER KRRQLEGTSE LPVKKAFRRG DWEKEREKKY LQEKQQKDEQ<br>RELKRRKLEE ERLKYEEKKL RISRLANKES SRNEELLTET TTPSPAVKAS<br>PASTKLSVSE NDRLSIPEIT KDNLTLETEI AKLREMKEPI RLFGESEET<br>IQRYYSLLKY KKLEEIENEL LTKGVETIDF EHATTTKPKV SKQVVAFLQH<br>GIRIWDNFLS SKSINSFESS ESQMQLKIFR QAKQDLVLI QLIVDEALND<br>DIFKSIAEIC YRCQKHEFVK ANDMYLRLTI GNAPWPIGVT MVGIHERSAH<br>QRLQANPSSN ILKDEKRRKC LQALKRFITF QERESSNLPE YTD   |
| <b>Source</b>            | Mammalian cell  |
| <b>Target Names</b>      | prp18   |
| <b>Protein Names</b>     | Recommended name: Pre-mRNA-splicing factor 18   |
| <b>Expression Region</b> | 1-343   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |