



# Recombinant Schizosaccharomyces pombe Enhancer of mRNA-decapping protein 3 (edc3)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-BP007393SXV   |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.   |
| <b>Uniprot No.</b>       | O94752  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)   |
| <b>Purity</b>            | ≥85% (SDS-PAGE)   |
| <b>Sequence</b>          | MSVADFYGSN VEVLLNNSDK ARGVITNFDS SNSILQLRLA NDSTKSIVTK<br>DIKDLRILPK NEIMPKNGTK SPSTNSTKLNK SAETYSSKNK WSMDCDEEFD<br>FAANLEKFDK KQVFAEFREK DKKDPAKLLV SHNKSPNRNY HHKQNVLGPS<br>VKDEFVDLPS AGSQINGIDA VLSSSSNGHV TPGSKKGSRE TLKKKPFVDE<br>NIPAEHHTTT GDILKPITPE QLSQGIALAI AKTSTDIVVE NAAQLLSQFV<br>FSVLGGHKRL SSRNHNSQPL VCILVGSHDH ASAAVAAGRR LCAIGIKVVL<br>RLLTPFNVDN RQLLMFQAAG GYIPTENFDQ FLNKLTSPIE LVVDVLTGFH<br>PSIDKNSHAL IQWANDLNVL ILSVDIPSGY TVQKKNTAIL PKWTLALGAV<br>TTTLAQAALV KQAAGVSVFV GNLGTGSQTW AELGILESQV TGQYLAQISC<br>TSTN |
| <b>Source</b>            | Baculovirus   |
| <b>Target Names</b>      | edc3  |
| <b>Protein Names</b>     | Recommended name: Enhancer of mRNA-decapping protein 3  |
| <b>Expression Region</b> | 1-454   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |