



Recombinant Schizosaccharomyces pombe Enhancer of mRNA-decapping protein 3 (edc3)

Product Code	CSB-MP007393SXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O94752
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MSVADFYGSN VEVLLNNSDK ARGVITNFDS SNSILQLRLA NDSTKSIVTK DIKDLRILPK NEIMPKNGTK SPSTNSTKLNK SAETYSSKNK WSMDCDEEFD FAANLEKFDK KQVFAEFREK DKKDPAKLLV SHNKSPNRNY HHKQNVLGPS VKDEFVDLPS AGSQINGIDA VLSSSSNGHV TPGSKKGSRE TLKKKPFVDE NIPAEHHTTT GDILKPITPE QLSQGIALAI AKTSTDIVVE NAAQLLSQFV FSVLGGHKRL SSRNHNSQPL VCILVGSHDH ASAAVAAGRR LCAIGIKVVL RLLTPFNVDN RQLLMFQAAG GYIPTENFDQ FLNKLTSPIE LVVDVLTGFH PSIDKNSHAL IQWANDLNVL ILSVDIPSGY TVQKKNTAIL PKWTLALGAV TTTLAQAALV KQAAGVSVFV GNLGTGSQTW AELGILESQV TGQYLAQISC TSTN
Source	Mammalian cell
Target Names	edc3
Protein Names	Recommended name: Enhancer of mRNA-decapping protein 3
Expression Region	1-454
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.