



# Recombinant Escherichia coli Phosphoserine aminotransferase (serC)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-MP018838ENT                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Uniprot No.</b>       | B1X847                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Immunogen Species</b> | Escherichia coli (strain K12 / DH10B)                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Purity</b>            | >85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Sequence</b>          | MAQIFNFSSG PAMPLAEVLK QAQQELRDWN GLGTSVMEVS HRGKEFIQVA<br>EEAEKDFRDL LNVPSNYKVL FCHGGGRGQF AAVPLNILGD KTTADYVDAG<br>YWAASAIKEA KKYCTPNVFD AKVTVDGLRA VKPMREWQLS DNAAYMHYCP<br>NETIDGIAID ETPDFGADV VVAADFSSSTIL SRPIDVSRYG VIYAGAQNKI<br>GPAGLTIVIV REDLLGKANI ACPSILDYSI LNDNGSMFNT PPTFAWYLSG<br>LVFKWLKANG GVAEMDKINQ QKAELLYGVI DNSDFYRNDV AKANRSRMNV<br>PFQLADSALD KLFLEESFAA GLHALKGHRV VGGMRASIYN AMPLGKAL<br>TDFMVEFERR HG |
| <b>Source</b>            | Mammalian cell                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Target Names</b>      | serC                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Protein Names</b>     | Recommended name: Phosphoserine aminotransferase EC= 2.6.1.52<br>Alternative name(s): Phosphohydroxythreonine aminotransferase Short name=<br>PSAT                                                                                                                                                                                                                                                                                 |
| <b>Expression Region</b> | 1-362                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                        |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                 |