



# Recombinant Callicebus moloch Serpin B5 (SERPINB5)

<b>Product Code</b>	CSB-YP021073CYE
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	B1MTB7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Plecturocebus moloch (Dusky titi monkey) (Callicebus moloch)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDALQLANSA FAVDMFKQLC EKEPVGNVLF SPICLSTSLA LAQVGAKGDT ANEIGQVLHF ENVKDVPFGF QVTSDVNKL SSFYSLKLIK RLYVDKSLNL STEFISSTKR PYAKELETVD FKDKLEETKG QINNSIKDLT DGHFENILAD NSVSNQTKIL VVNAAYFVGK WMKKFPESET KECPRVNKT DTKPVQMMNI EATFCMGNID SIDCKIIELP FQNKHLSMFI LLPKDVEDES TGLEKIEKQL NSEALAQWTN PSTMANAKVK LSIPKFKVEK IIDPKASLEN LGLKRIFSED TSDFGMSET KGVALSNVIH KVCLEITEDG GDSIEVPGAR ILQHKDELNA DHPFVYIIRH NKTRNIIFFG KFCSP
<b>Source</b>	Yeast
<b>Target Names</b>	SERPINB5
<b>Protein Names</b>	Recommended name: Serpin B5 Alternative name(s): Maspin Peptidase inhibitor 5 Short name= PI-5
<b>Expression Region</b>	1-375
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.